

Profiling epitope-specific antibody responses against malaria utilizing high-density peptide arrays

Kirsten Heiss^{1*}, Thomas Jaenisch^{2*}, Nico Fischer², Carolin Geiger², Gregor Jainta¹, Leszek Rychlewski³, Ali Sie⁴, Boubacar Coulibaly⁴, Lucjan S. Wyrwicz⁵, Frank Breittling⁶, and Felix F. Loeffler⁷






¹PEPPERPRINT GmbH, Heidelberg, Germany; ²Center for Infectious Diseases & Heidelberg Institute of Global Health, Heidelberg University Hospital, Germany; ³BioInfoBank Institute, Poznań, Poland; ⁴Centre de Recherche en Santé de Nouna, Nouna, Burkina Faso; ⁵Department of Oncology and Radiotherapy, M. Skłodowska Curie Memorial Cancer Center, Warsaw, Poland; ⁶Institute of Microstructure Technology, Karlsruhe Institute of Technology, Karlsruhe, Germany; ⁷Max Planck Institute of Colloids and Interfaces, Potsdam, Germany

PEPPERPRINT
A NEW DIVERSITY

Abstract

Antibody-mediated defense mechanisms play an important role in combating malaria infection attacking the parasite at multiple life cycle stages. By reducing parasitemia and clinical symptoms, humoral responses are the key immune effector mechanisms at the pathogenic asexual blood stage. Intensive research focusses on the identification of target antigens of protective anti-plasmodial immune responses, but the entire picture is still incomplete.

High-density peptide microarrays represent an emerging tool, which can be applied to screen a large diversity of linear and conformational antibody epitopes. The technology provides a multifaceted application spectrum, where antibody reactivities to tens of thousands of peptides can be monitored simultaneously in a single assay. We used the peptide microarray technology to better characterize epitopes associated with naturally acquired immunity (NAI). We analyzed the peptide array profiles of serum samples from individuals living in a malaria-endemic area and compared the responses to malaria-naïve individuals. The results show distinct antibody patterns according to immune status for peptides predicted by applying different algorithms, as well as for peptides derived from well-known vaccine candidates. We demonstrated that peptide arrays can be applied as a novel screening method for the identification of new immunogenic antibody epitopes. The discovery of the antigens and/or antigenic motifs that are responsible for protective immune responses will ultimately create new opportunities for vaccine development.



Summary & Conclusion

- First study using high-density peptide arrays to profile anti-plasmodial antibody responses in NAI against linear epitopes
- Identification of novel, highly immunogenic epitopes in both well-known vaccine candidates and previously uncharacterized antigens
- The peptide-based approach allows to unravel differential antibody binding to specific epitopes which might be missed when using the entire antigen