

PCR Plates & Accessories

Features

- USP class VI certified polypropylene plates
- Thin-walled for optimal heat transfer
- Free of detectable DNase, RNase, DNA, PCR Inhibitors, ATP, Endotoxins
- Designed to minimize sample loss
- Broad compatibility with RT-PCR instruments

i-labware
more than just equipment



Ordering Information

Cat No.	Description	Package Quantity
C20010	0.2ml PCR 96 well plate, clear, non-skirted	10/bag, 5 bags/case
C20020	0.1ml qPCR 96 well plate, clear, sub-semi skirted	10/bag, 5 bags/case
C20021	0.1ml qPCR 96 well plate, low profile, white, semi skirted	10/bag, 5 bags/case
C20022	0.2ml qPCR 96 well plate, clear, sub-semi skirted	10/bag, 5 bags/case
C30010	Storage sealing membrane	10/bag, 10 bags/case
C30011	Aluminum sealing membrane	10/bag, 10 bags/case
C30020	PCR 96 well plate sealing membrane	10/bag, 10 bags/case
C30021	qPCR 96 well plate optical sealing membrane	10/bag, 10 bags/case

i-Labware Pte Ltd

102E Pasir Panjang Road, #08-03 Citilink Warehouse Complex, Singapore 118529

enquiry@i-labware.com | www.i-labware.com

Instrument Compatibility for PCR Plate

	Model	C20010	C20020	C20021	C20022
	Standard Well (0.2ml, height 21mm)	✓			✓
	Low Profile Well (0.1ml, height 16mm)		✓	✓	
Applied Biosystems	7000				✓
	7300				✓
	7500				✓
	7900HT Fast (with tray) (0.1ml type)		✓		
	7900HT (with tray) (0.2ml type)				✓
	StepOne Plus		✓		
	7500Fast		✓		
	ViiA7 Fast (0.1ml model)		✓		
	ViiA7 (0.2ml model)				✓
	QuantStudio Fast 96-well Block (0.1ml model)		✓		
	QuantStudio 96-well Block (0.2ml model)				✓
	5700				✓
	7700				✓
	310/3100/3130				✓
	3700/3730				✓
Roche	LightCycler 480			✓	
	LightCycler 480 II			✓	
	LightCycler 96			✓	
Bio-Rad	iCycler iQ	✓			
	iQ5	✓			
	MyiQ	✓			
	Chromo4	✓			
	Opticon I			✓	
	Opticon II			✓	
	CFX96			✓	
Others	Stratagene Mx4000/3000P/3005P	✓			
	Eppendorf MasterCycle ep Realplex	✓			
	Biometra Toptical Thermocycler	✓			