


Pipette Tips

Features

- Unique conical design, compatible with most pipette brands (Biohit, Eppendorf, Gilson, Rainin, Finnpipette)
- RNase, DNase and Pyrogen Free
- Available in multiple formats (bags, racks, refills, sterile, non-sterile, filtered, non-filtered, low retention)
- High density polyethylene filters
- Tips and racks are fully autoclavable
- Extended length tips greatly reduces the possibility of sample-pipettor contact
- Low Retention tips ensure near complete sample recovery

i-labware
more than just equipment



i-Labware Pte Ltd

102E Pasir Panjang Road, #08-03 Citilink Warehouse Complex, Singapore 118529
enquiry@i-labware.com | www.i-labware.com

Ordering Information

Universal Fit Pipette Tips

Cat No.	Description	Package Quantity
C41011	0.1-10 ul, clear tips, bulk pack	1000/bag, 10 bags/case
C41031	0.1-10 ul, clear tips, sterile, rack pack	96/rack, 50 racks/case
C41051	0.1-10 ul, clear tips, sterile, refill pack	480/rack, 10 racks/case
C41212	1-200 ul, yellow tips, bulk pack	1000/bag, 10 bags/case
C41232	1-200 ul, yellow tips, sterile, rack pack	96/rack, 50 racks/case
C41252	1-200 ul, yellow tips, sterile, refill pack	480/rack, 10 racks/case
C41313	100-1000 ul, blue long tips, bulk pack	1000/bag, 10 bags/case
C41333	100-1000 ul, blue long tips, sterile, rack pack	96/rack, 50 racks/case
C41353	100-1000 ul, blue long tips, sterile, refill pack	480/rack, 10 racks/case

Universal Fit Filter Pipette Tips

Cat No.	Description	Package Quantity
C45031	0.1-10 ul, clear filter tips, sterile, rack pack	96/rack, 50 racks/case
C45532	1-200 ul, yellow filter tips, sterile, rack pack	96/rack, 50 racks/case
C45333	100-1000 ul, blue filter tips, sterile, rack pack	96/rack, 50 racks/case

Low Retention - Universal Fit Pipette Tips

Cat No.	Description	Package Quantity
C41011-L	0.1-10 ul, clear tips, low retention, bulk pack	1000/bag, 10bags/case
C41031-L	0.1-10 ul, clear tips, low retention, sterile, rack pack	96/rack, 50 racks/case
C41212-L	1-200 ul, yellow tips, low retention, bulk pack	1000/bag, 10bags/case
C41232-L	1-200 ul, yellow tips, low retention, sterile, rack pack	96/rack, 50 racks/case
C41313-L	100-1000 ul, blue long tips, low retention, bulk pack	1000/bag, 10bags/case
C41333-L	100-1000 ul, blue long tips, low retention, sterile, rack pack	96/rack, 50 racks/case

Low Retention - Universal Fit Filter Pipette Tips

Cat No.	Description	Package Quantity
C45031-L	0.1-10 ul, clear filter tips, low retention, sterile, rack pack	96/rack, 50 racks/case
C45532-L	1-200 ul, yellow filter tips, low retention, sterile, rack pack	96/rack, 50 racks/case
C45333-L	100-1000 ul, blue long filter tips, low retention, sterile, rack pack	96/rack, 50 racks/case